

Product Datasheet

Staphylococcus aureus Alpha-hemolysin antibody (orb243684)

Description	Staphylococcus aureus Alpha-hemolysin antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Staphylococcus aureus Alpha-hemolysin protein (27-319AA)
Target	hly
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q2G1X0
Expiration Date	12 months from date of receipt.