

Product Datasheet

Staphylococcus aureus Staphopain B antibody (orb243688)

Description	Staphylococcus aureus Staphopain B antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Staphylococcus aureus Staphopain B protein (220-393AA)
Target	sspB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99V46
Expiration Date	12 months from date of receipt.