

Product Datasheet

Sus scrofa Pulmonary surfactant-associated protein B antibody (orb243696)

Description	Sus scrofa Pulmonary surfactant-associated protein B antibody
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Pig Pulmonary surfactant-associated protein B protein (1-79AA)
Target	SFTPB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15782
Expiration Date	12 months from date of receipt.