

Product Datasheet

CCDC88B antibody (orb2437)

Description	CCDC88B antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	KLH conjugated synthetic peptide derived from human CCDC88B/BRLZ/HkRP3 (1271-1476/1476 aa)
Target	CCDC88B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	79/165 kDa
Uniprot ID	A6NC98
Expiration Date	12 months from date of receipt.