

Product Datasheet

Albumin antibody (orb243719)

Description	Albumin antibody
Species/Host	Guinea pig
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Serum Albumin purified from Rabbit serum (Native Protein)
Target	IgG
Preservatives	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49065
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.