

Product Datasheet

Albumin antibody (Biotin) (orb243722)

Description	Albumin antibody (Biotin)
Species/Host	Guinea pig
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Serum Albumin purified from Rabbit serum (Native Protein)
Target	IgG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49065
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.