

Product Datasheet

Prealbumin protein antibody (orb243735)

Description	Prealbumin protein antibody
Species/Host	Guinea pig
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human Prealbumin protein (Native Protein)
Target	TTR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02766
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.