

Product Datasheet

Ovotransferrin antibody (HRP) (orb243740)

Description	Ovotransferrin antibody (HRP)
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Chicken Ovotransferrin (Native Protein)
Target	Ovotransferrin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02789
Dilution Range	ELISA: 1:200
Expiration Date	12 months from date of receipt.