

Product Datasheet

Chorionic somatomammotropin hormone antibody (orb243743)

Description Chorionic somatomammotropin hormone antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen Recombinant Human Chorionic somatomammotropin hormone 1 protein (27-217AA)

Target CSH1

Preservatives Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only

Isotype IgG

Clonality Polyclonal

Uniprot ID [P0DML2](#)

Expiration Date 12 months from date of receipt.



Western blot analysis of A549 (Lane 1), u...