



Product Datasheet

Inhibin beta B chain antibody (orb243758)


Description	Inhibin beta B chain antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Inhibin beta B chain protein (295-405AA)
Target	INHBB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09529
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200, IF/ICC:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Immunohistochemical staining of human liver tissue.



Immunofluorescence analysis of PC-3 cells.