

Product Datasheet

Mycobacterium tuberculosis 85-A antibody (FITC) (orb243768)

Description	Mycobacterium tuberculosis 85-A antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Mycobacterium tuberculosis Diacylglycerol acyltransferase/mycolyltransferase Ag85A protein (53-331AA)
Target	fbpA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WQP3
Expiration Date	12 months from date of receipt.