

Product Datasheet

IgG Fc fragment antibody (FITC) (orb243780)

Description	IgG Fc fragment antibody (FITC)
Species/Host	Rabbit
Reactivity	Goat
Conjugation	FITC
Immunogen	Goat IgG Fc fragment
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.