

Product Datasheet

Mycobacterium tuberculosis ESAT-6-like protein EsxB antibody (Biotin) (orb243786)

Description	Mycobacterium tuberculosis ESAT-6-like protein EsxB antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mycobacterium tuberculosis ESAT-6-like protein EsxB protein (2-100AA)
Target	esxB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WNK4
Expiration Date	12 months from date of receipt.