



Product Datasheet

Pro-opiomelanocortin antibody (orb243791)

Description	Pro-opiomelanocortin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Peptide sequence from Human Pro-opiomelanocortin protein
Target	POMC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01189
Dilution Range	WB:1:1000-5000, IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 whole ce...



Immunohistochemical staining of human br...