

Product Datasheet

IgG heavy chain antibody (Biotin) (orb243831)

Description	IgG heavy chain antibody (Biotin)
Species/Host	Guinea pig
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Human IgG heavy chain (Native Protein)
Target	IgG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.