



Product Datasheet

Inosine-5'-monophosphate dehydrogenase 2 antibody (orb243836)


Description	Inosine-5'-monophosphate dehydrogenase 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Inosine-5'-monophosphate dehydrogenase 2 protein (5-514AA)
Target	IMPDH2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12268
Dilution Range	WB:1:500-2000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of Hela cells...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...