

Product Datasheet

Inosine-5'-monophosphate dehydrogenase 2 antibody (HRP) (orb243837)

Description	Inosine-5'-monophosphate dehydrogenase 2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Inosine-5'-monophosphate dehydrogenase 2 protein (5-514AA)
Target	IMPDH2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12268
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.