

Product Datasheet

Mutant hemoglobin alpha 1 globin chain antibody (Biotin) (orb243841)

Description	Mutant hemoglobin alpha 1 globin chain antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Hemoglobin subunit alpha protein (2-142AA)
Target	HBA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P69905
Expiration Date	12 months from date of receipt.