

Product Datasheet

Vasopressin-neurophysin 2-copeptin antibody (HRP) (orb243847)

Description	Vasopressin-neurophysin 2-copeptin antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Vasopressin-neurophysin 2-copeptin protein (126-164AA)
Target	AVP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01185
Expiration Date	12 months from date of receipt.