

Product Datasheet

TNF alpha antibody (orb243864)

Description	TNF alpha antibody
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Bovine TNF alpha protein
Target	TNF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q06599
Expiration Date	12 months from date of receipt.