



Product Datasheet

RBP4 antibody (orb243868)

Description	RBP4 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Retinol-binding protein 4 protein
Target	RBP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P02753
Dilution Range	Recommended dilution:WB:1:500-1:5000,IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Human positive ...



Western blot analysis of RBP4 transfecte...