


Product Datasheet

IL6 antibody (orb243911)

Description	IL6 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Interleukin-6
Target	IL6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P05231
Dilution Range	Recommended dilution:WB:1:500-1:5000,IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pr...