

Product Datasheet

Folic acid antibody (orb243913)

Description	Folic acid antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Folic acid-BSA conjugate
Target	Folic acid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	ELISA: 0.5-40.5 ppb
Expiration Date	12 months from date of receipt.