

Product Datasheet

The metabolite of furazolidone antibody (orb243915)

Description	The metabolite of furazolidone antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	The metabolite of furazolidone-BSA conjugate
Target	The metabolite of furazolidone
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.