

Product Datasheet

TFF3 antibody (orb243923)

Description	TFF3 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Trefoil factor 3
Target	TFF3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q07654
Dilution Range	IHC-P:1:50-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.