

Product Datasheet

Syndecan 1 antibody (orb243927)

Description	Syndecan 1 antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Rat Syndecan-1 protein
Target	Sdc1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P26260
Expiration Date	12 months from date of receipt.