

Product Datasheet

IL1F10 antibody (orb243931)

Description	IL1F10 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human IL1F10 protein
Target	IL1F10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q8WWZ1
Dilution Range	Recommended dilution:IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.