

Product Datasheet

GST antibody (orb243938)

Description	GST antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant GST Protein
Target	GST
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	Recommended dilution:WB:1:500-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of using GST antib...