



Product Datasheet

Ki-67 antibody (PE) (orb243941)

Description Mouse monoclonal antibody against Ki-67 conjugated to


Flow
cytometry
analysis of
human
periphe...

Reactivity	Bovine, Human
Conjugation	PE
Tested Applications	FC
Immunogen	Nuclei of the Hodgkin lymphoma cell line L428
Target	Ki-67
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P46013
Entrez	4288
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.
Expiration Date	12 months from date of receipt.