

Product Datasheet

Human ASPH protein (orb244034)

Description	Recombinant human Aspartyl/asparaginyl beta-hydroxylase
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Partial of the full length of Isoform 6 of His-SUMO-tag and expression region is 75-270aa
Protein Sequence	FDLVDYEEVLGKLGIDADGDGDFDVEDDAKVLLGLKERSTSEPAVPPEEAEPHTE PEEQVPVEAEPQNIIDEAKEIQSLLHEMVHAEHETEHSYHVEETVSQDCNQDM EEMMSEQENPDSSEPVEDERLHHDTDDVTYQVYEEQAVYEPLENEGIEITEVT APPEDNPVEDSQVIVEEVSIFPVEEQQEVPPDT
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	38.2 kDa
Uniprot ID	Q12797
Expiration Date	6 months from date of receipt.