

Product Datasheet

Human EYA2 protein (orb244192)

Description	Recombinant human Eyes absent homolog 2
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Isoform 2 of His-SUMO-tag and expression region is 1-514aa
Protein Sequence	MVELVISPLTVNSDCLDKLFNRADAAVWTLSDRQGITSAPLRVSQFLFSRSCP RVLPRQPSTAMAAYGQTQYSAGIQATPYTAYPPPAQAYGIPSFSTSGTQSPYT YQMHGTTGFYQGGNGLGNAAGFGSVHQDYPSYPGFPQSQYPQYYGSSYNPPY VPASSICPSPLSTSTYVLQEASHNVPNQSSESLAGEYNTHNGPSTPAKEGDTDRP HRASDGKLRGRSKRSSDPSPAGDNEIERVFVWDLDETIIIFHSLLTGTFFASRYGK DTTTSVRIGLMMEEMIFNLADTHLFFNDLEDQDIHVDDVSSDDNGQDLSTYNF SADGFHSSAPGANLCLGSGVHGGVDWMRKLAFRYRRVKEMYNTYKNNVGGGLI GTPKRETWLQLRAELEALTDLWLTHSLKALNLINSRPNVNVLTQTLPALAKV LLYGLGSVFPIENIYSATKTGKESCFERIMQRFGRKAVYVVGIDGVEEEQGAKKH NMPFWRISCHADLEALRHALELEYL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	72.7 kDa
Uniprot ID	O00167
Expiration Date	6 months from date of receipt.