

Product Datasheet

Human GLCE protein (orb244232)

Description	Recombinant human D-glucuronyl C5-epimerase
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Expression region is 29-617aa and His-tag protein
Protein Sequence	NKCSSDKAIQFPRSSSGFRVDGFEKRAAASESNMNHVAKQQSEEAFFPQEQ QKAPPVVGGFNSNVGSKVLGLKYEEIDCLINDEHTIKGRREGNEVFLPFTWVEKY FDVYGKVVQYDGYDRFEFSHSYSKVYAQRAPYHPDGVFMSFEGYNVEVRDRVK CISGVEGVPLSTQWGPQGYFYPIQIAQYGLSHYSKNLTEKPPHIEVYETAEDRDKN KPNDWTVPKGCFMANVADKSRFTNVKQFIAPETSEGVSLQLGNTKDFIISFDLKF LTNGSVSVVLETTTEKNQLFTIHVYSNAQLIAFKERDIYYGIGPRTSWSTVTRDLVT DLRKGVGLSNTKAVKPTKIMPKKVRLIAKGKGLDNITISTTAHMAAFFAASDW LVRNQDEKGGWPIMVTRKLGEFGKSLEPGWYSAMAQQAISTLVRAYLLTKDHI FLNSALRATAPYKFLSEQHGVKAVFMNKHWDYEEYPTTPSSFVLNGFMYSLIGLY DLKETAGEKLGKEARSLYERGMESLKAMLPLYDTGSGTIYDLRHFMLGIAPNLAR WDYHTTHINQLQLLSTIDESPVFKEFVKRWKSYLKGSRAKHN
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	70.9 kDa
Uniprot ID	094923
Expiration Date	6 months from date of receipt.