

## Product Datasheet

### Human HJURP protein (orb244271)

**Description**

Recombinant human Holliday junction recognition protein

**Endotoxins**

Not test.

**Preservatives**

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

**Form/Appearance**

Liquid or Lyophilized powder

**Storage**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

**Tag**

N-terminal 6xHis-SUMO-tagged

**Note**

For research use only

**Application notes**

Full length of His-SUMO-tag and expression region is 1-748aa

**Protein Sequence**

MLGTLRAMEGEDVEDDQLLQKLASRRRFQRRMQRLIEKYNQPFEDTPVVQMA  
TLTYETPQGLRIWGGRLIKERNEGEIQDSSMKPADRTDGSVQAAAWGPELPSHR  
TVLGADSKSGEVDATSDQEESVAWALAPAVPQSPLKNELRRKYLTVQVDILLQGA  
EYFECAGNRAGRDRVVTPLPSLASPAVPAPGYCSRISRKSPGDPKAPASSPREW  
DPLHPSSTDMALVPRNDSLQETSSSSFLSSQPFEDDDICNVITISDLYAGMLHS  
MSRLLSTKPSSIISTKTFIMQNWNSRRRHRYKSRMNKTYCKGARRSQRSSKENFI  
PCSEPVKGTGALRDCKNVLDVSCRKTGLKLEKAFLEVNRPQIHKLDPWKERKV  
TPSKYSSLIYFDSSATYNLDEENRFRTLKWLISPVKIVSRPTIRQGHGENRQREIEI  
RFDQLHREYCLSPRNQPRRMCLPDSWAMNMYRGGPASPGGLQGLETRRLSLPS  
SKAKAKSLSEAFENLGKRSLEAGRCLPKSDSSSLPKTNPTHSATRPQQTSDLHV  
QGNSSGIFRKSVSPSKTSLVDPKEVPGHGRNRYDEIKEEFDKLHQKYCLKSPGQ  
MTVPLCIGVSTDKASMEVRYQTEGFLGKLNPDPHFQGFQKLPSSPLGCRKSLLG  
STAIEAPSSTCVARAITRDGTRDHFPAKRPRLSEPGSGRQGNLSLGASDGVDN  
TVRPGDQGSSSQPNSEERGENTSYRMEEKSDFMLEKLETKSV

**Purity**

Greater than 90% as determined by SDS-PAGE.

**Source**

E.coli

**MW**

99.5 kDa

**Uniprot ID**
[Q8NCD3](#)
**Expiration Date**

6 months from date of receipt.