

Product Datasheet

Rat Hpx protein (orb244277)

Description	Recombinant rat Hemopexin
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full length of His-tag and expression region is 24-460aa
Protein Sequence	NPLPAAHETVAKGENGTKPDSVDIEHCSDAWSFDATTMDHNGTMLFFKGEFVW RGHSGIRELISERWKNPVTSDAAFRGPDSVFLIKEDKVVVYPPEKKENGYPKLF QEEFPGIPYPPDAAVECHRGECQSEGVLFQGNRKWFWDFA TRTQKERSWPAV GNCTAALRWLERYYCFOGNKFLRFNPVTGEVPPRYPLDARDYFISCPGRGHGKL RNGTAHGNSTHPMHSRCNADPGLSALLSDHRGATYAFSGSHYWRLDSSRDGW HSWPIAHHWPQGPSAVDAAFSWDEKVYLIQGTQVYVFLTKGGNNLVSGYPKRL EKELGSPPGISLDTIDAAFSCPGSSKLYVTSGRRLWWLDLKSGAQATWAELSWP HEKVDGALCLEKSLGPYSCSSNGPNLFFIHGPNLYCYSSIDKLNAAKSLPQPQKV NSILGCSQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	52.9 kDa
Uniprot ID	P20059
Expiration Date	6 months from date of receipt.