

Product Datasheet

Human KRT4 protein (orb244333)

Description	Recombinant human Keratin, type II cytoskeletal 4
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Full length of His-SUMO-tag and expression region is 1-534aa
Protein Sequence	MIARQQCVRRGGPRGFSCGSAIVGGGKRGAFSSVSMSSGGAGRCSSGGFGSRSLY NLRGNKSISMSVAGSRQACFGGAGGFGTGGFGAGGFGAGFGTGGFGGGFG GSFSGKGGPGFPVCPAGGIQEVNTINQSLLTPLHVEIDPEIQKVRTEEREQIKLLNN KFASFIDKVQFLEQQNKVLETKWNLLQQQTTTTSSKNLEPLFETYLSVLRKQLDT LGNDKGRLQSELKTMQDSVEDFKTYEEEINKRTAAENDFVVLKKDVEDAAYLNK VELEAKVDLSLNDEINFLKVLYDAELSQMQTHVSDTSVVLMDNNRNLDLDSIAE VRAQYEEIAQRSKAEAEALYQTKVQQLQISVDQHGDNLKNKSEIAELNRMIRQL RAEINIKKQCQTLQVSVADAEQRGENALKDAHSKRVELEAALQQAKEELARML REYQELMSVKLALDIEIATYRKLLGEEYRMSGECQSAVSISVSVSGSTSTGGISGG LGSGSGFGLSSGFGSGSGSGFGFGGSGVSGSSSSKIISTTTLNKRR
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	73.3 kDa
Uniprot ID	P19013
Expiration Date	6 months from date of receipt.