

Product Datasheet

Human LRBA protein (orb244354)

Description	Recombinant human Lipopolysaccharide-responsive and beige-like anchor protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Partial of the full length of 2-2863aa of His-tag and expression region is 1-251aa
Protein Sequence	MASEDNRVSPSPPTGDDGGGGGREETPTEGGALSLKPGLPPIRGIRMKFAVLTGL VEVGEVSNRDIVETVFNLLVGGQFDLEMNFIIQEGESINCMVDLLEKCDITCQAE VWSMFTAILKKSIRNLQVCTEVGLVEKVLGKIEKVDNMIADLLVDMGLVLSYNL TVRELKLFSSKLQGDKGRWPPHAGKLLSVLKHMPQKYGPDAFFNFPKSAIAA LPPIAKWPYQNGFTFHTWLRMDPVNNINVDKDK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	31.6 kDa
Uniprot ID	P50851
Expiration Date	6 months from date of receipt.