

Product Datasheet

Human NXF3 protein (orb244418)

Description	Recombinant human Nuclear RNA export factor 3
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Full length of His-SUMO-tag and expression region is 1-531aa
Protein Sequence	MSLP SGHTTGHTDQVVQRRARCWDIYQRRFSSRSEPVNPGMHSSSHQQQDG DAAMHGAHMDSPVRYTPYTISPYNRKGSFRKQDQTHVNMEREQKPPERRMEG NMPDGTLGSWFKITVPFGIKYNEKWLLNLIQNECSVPFVPVEFHYENMHASFFVE NASIAYALKNVSGKIWDEDNEKISIFVNPAGIPHFVHRELKSEKVEQIKLAMNQOC DVSQEALDIQRLPFYPMVNRTKMASNPRKCMAASLDVHEENIPTVMSAGEM DKWKGIEPGEKCADRSPVCTTFSDTSSNINSILELFPKLLCLDGQQSPRATLCGTE AHKRLPTCKGSFFGSEMLKNLVLQFLQYYLIYDSGDRQGLLSAYHDEACFSLSI PFNPEDSAPSSFCFFKDSRNILKDPYLRGELLKHTKLDIVDSLSALPKTQHDL SFLVDMWYQTEWMLCFVNGVFKEVEGQSQGSVLAFTRTFIATPGSSSSLCIVN DKLFVRDTS HQGTQSALFTLVPTAFSSSVPAFSQEQQKMLPS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	76.1 kDa
Uniprot ID	Q9H4D5
Expiration Date	6 months from date of receipt.