

Product Datasheet

Mouse Pitx3 protein (orb244442)

Description	Recombinant mouse Pituitary homeobox 3
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	Full length of GST-tag and expression region is 1-302aa
Protein Sequence	MEFGLLGAEARSPALSLSDAGTPHPPLPEHGCKGQEHSDSEKASASLPGGSPEDGSLKKKQRRQRTHFTSQQLQELEATFQRNRYPD MSTREEIAVWTNLTEARVRVWFKNRRRAKWRKRERSQQAELCKGGFAAPLGGLVPPYEEVYPGYSYGNWPPKALAPPLAAKTFFAFNSVNVGPLASQPVFSPSSIAASMVPSAAAAAPGTVPGPGALQGLGGAPPGLAPAAVSSGAVSCPYASAAAAAASSPYVYRDPCNSSLASLRLKAKQHASFSPAVPGPPPAANLSPCQYAVERPV
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	58.7 kDa
Uniprot ID	O35160
Expiration Date	6 months from date of receipt.