

Product Datasheet

Mouse PRKDC protein (orb244459)

Description	Recombinant mouse DNA-dependent protein kinase catalytic subunit
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	PI3K/PI4K domain of His-tag and expression region is 3747-4015aa
Protein Sequence	EYPFLVKGGEDLRQDQRIEQIFEVMNAILSQDAACSQRNMQLRITYRVVPMTSRL GLIEWIENTMTLKDLLLSNMSQEEKVANNSDPKAPIRDYKDWLMKVSGKSDAGA YVLMYSRANRTETVVAFRRRRESQVPPDLLKRAVVKMSTSPEAFLALRSHFASSHA LLCISHWLLGIGDRHLNFMVAMETGSVIGIDFGHAFGSATQFLPVPELMPFRLT RQFVSLMLPMKETGLMCTVMVHALRAFRSCAGLLTDTMEIFVKEPSFDWKS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	34.6 kDa
Uniprot ID	P97313
Expiration Date	6 months from date of receipt.