

Product Datasheet

Chicken VTG1 protein (orb244652)

Description	Recombinant chicken Vitellogenin-1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Vitellogenin of His-tag and expression region is 24-633aa
Protein Sequence	FGENKVYTYNYESILFSGIPEKGLARTGIRIRSEVEISGIGPKLCLIRIHSIEAAEYNGI WPTSSFSRSLKLTQALTGQLSIPIKFEYSNGHVGNLMAPDSVSDDGILNIYRGILNI LELSLKKMQHSYSIQEAGIGGICNTTYAIQENKKANLVDVTKSKDLNSCEEKVQV VTGSAYTQPCQTCQQRNKNRATATYNYKIKYTHNEAVITQAEVEEVHQFTPFHE ITGGNAIVEARQKLALIEVQKQVAEVPPEFQKRGSLQYQFGSELLQLPVHLFKIK DVERQIEERLQDLVETTYEQLPSDAPAKALKLMHLLRAANEENYESVWKQFSSR PAYRRYLLDLLPAAASHRSLRFLRHKMERQELTNWEIAQTVLVALHSSSPTQEV EATLIVKKHCPSSSVLRKVCLLSYASLCHKRCSSPYSCSECLQVFHVFAGEALG KSNIEEVLLALKALGNVGHPPASIKHIKKFLPGYAAGASELPLKVHETAVMALKSIG MRDPQMVQAITLEIFLNHKKIHPRIIRMLAAVVLETKPGLPILMILVDAVLKEPSMQV ASFIYSHLRALGRSTAPDLQMMASACRMAVRALSPKFDMSGYQFSKVFRFSMFK EFLMSGLAACYFVLNNAGSLIPTMAV
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	75.6 kDa
Uniprot ID	P87498
Expiration Date	6 months from date of receipt.