

Product Datasheet

E.coli IFM1 protein (orb244696)

Description

Recombinant *Saccharomyces cerevisiae* (strain ATCC 204508 / S288c)

Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	This is GST-tag protein
Protein Sequence	KGKRRNQISKKEKPLNFSIPNYISVNKLANLLNCRVERLIKDLTALGFENITTTYIL SKEYVELILQEYNFALPNLSTSTNLDNVYDELKSPVNPKLLTKRAPVVTIMGHVDH GKTTIIDYLRKSSVVAQEHGGITQHIGAFQITAPKSGKKITFLDTPGHA AFLKMRER GANITDIIVLVSVESLMPQTLEAIKHAKNSGNEMIIAITKIDRIPQPKEREKKIEKV INDLIVQGIPVEKIGGDVQVIPISAKTGENMDLLEESIVLLSEVMDIRAENSPKTI AE GWIIESQVKKQVGNVATVLVKKGTLQKGKILICGNTFCIKNLIDDKGIPILKATPS YATEVLGWKDVPHVGDEVIQVKSEAIKKFISKRQDLIEVQKNSSIVEKLNEERAL AKEQHLNKELEHENTVQEHEQNTGPKLINYIIKCDVSGSAEAVSESISSLGNDEV RCNVISSSVGIPTESDLKMAQITESTILCFNLGNLPSEVINNRAGIKIKQYNVIYKLIE DVTETLTENLKPIFEKKIVSTVDVRETFDFRLKKKIIRIAGCKVNNGVIKNSLVQVV RGPNEVDVFDGKISTLKH NKDDVAEVSKGHECGITFESGFEGFKPGDKILVYENV RVPRYL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	97.5 kDa
Uniprot ID	P25038
Expiration Date	6 months from date of receipt.