

Product Datasheet

Human RALGPS1 protein (orb244842)

Description	Recombinant human Ras-specific guanine nucleotide-releasing factor RalGPS1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	MYKRNLMASVLVTSATPQGSSSSDSLEGQSCDYASKSYDAVVFDVLKVTPEEF ASQITLMDIPVFKAIQPEELASCGWSKKEKHS LAPNVVAFTRRFNQVSFWVVREI LTAQTLKIRAEILSHFVKIAKKLLELNNLHSLMSVVSALQSAPIFRLTKTWALLNRK DKTTFEKLDYLM SKEDNYKRTREYIRSLKMVPSIPYLGIIYLLDIYIDSAYPASGSIM ENEQRSNQMNILRIIADLQVSCSYDHLTTLPHVQKYLKSVRYIEELQKFVEDDNY KLSLRIEPGSSSPRLVSSKEDLAAM
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	50.7 kDa
Uniprot ID	Q5JS13
Expiration Date	6 months from date of receipt.