

Product Datasheet

Human TBCEL protein (orb244851)

Description	Recombinant human Tubulin-specific chaperone cofactor E-like protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	MDQPSGRSFMQVLCEKYSPENFPYRRGPGMGVHVPATPQGSPMKDRLNLPSVL VLNSCGITCAGDEKEIAAFCAHVSELDLSDNKLEDWHEVSKIVSNVPQLEFLNLS SNPLNLSVLERTCAGSFSGVRKLVLNNSKASWETVHMILQELPDLEELFLCLNDY ETVSCPSICCHSLKLLHITDNNLQDWTEIRKLGVMFPSLDTLVANNHLNAIEEPD DSLARLFPNLRISLHKSGLSQSWEDIDKLNSFPKLEEVRLGIPLLQPYTTEERRKL VIARLPVSKLNGSVVTDGEREDSERFFIRYYVDVPQEEVPFRYHELITKYGKLEPL AEVDLRPQSSAKVEVHFNDQVEEMSIRLDQTVAEKKQLKTLVQLPTSNMLLYYF DHEAPFGPEEMKYSSRALHSFGIRDGDGKIYVESKTK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	64.2 kDa
Uniprot ID	Q5QJ74
Expiration Date	6 months from date of receipt.