

Product Datasheet

Human PLD5 protein (orb244892)

Description	Recombinant human Inactive phospholipase D5
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	MGEDDGLSEKNCQNKCRALVENIPEGLNYSNAPFHLSLFQGWMLLNMAK KSVDIVSSHWDNLNHTHPSACQGQRLFEKLLQLTSQNIKLVSDVTADSKVLEAL KLKGAEVTYMNMATAYNGRLQSSFWIVDKQHVYIGSAGLDWQSLGQMKELGVI FYNCSCVLVDLQRIKALYSSKFKSRVPQTWSKRLYGVDNEKKLQLNETKSQ AFVSNPKLFCPKNRSFDIDAIYSVIDDAKQYVYIAVMDYLPISSTSTKRTYWPDL AKIREALVLRVSRVRLLSFWKETDPLTFNFISLKAICTEIANCSLKVKFFDLEREN ACATKEQKNHTFPRLNRNKYMVTDGAAYIGNFDWVGNDFTQNAAGTGLVINQAD VRNNRSIIKQLKDVFERDWYSPYAKTLQPTKQPNCSSLFKLKPLSNKTATDDTGG KDPRNV
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	66.8 kDa
Uniprot ID	Q8N7P1
Expiration Date	6 months from date of receipt.