

Product Datasheet

Human TBX15 protein (orb244898)

Description	Recombinant human T-box transcription factor TBX15
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	MSSMEEIQVELQCADLWKRFDIGTEMIITKAGRRMFAMRVKITGLDPHQYYI AMDIVPVDNKRYYRYVYHSSKWMVAGNADSPVPPRVYIHPDSLASGDTWMRQV VSFDKLLTNNELDDQGHILHSMHKYQPRVHVIRKDFSSDLSPTKPVVPGDGVK TFNFPETVFTTVTAYQNQQITRLKIDRNPFAGFRDSGRNRTGLEAIMETYAFWR PPVRTLTFFEDFTTMQKQGGSTGTSPPTTSSTGTPSPSASSHLLSPSCSPPTFHLA PNTFNVGCRESQLCNLNSDYPPCARSNMAALQSYPLSDSGYNRLQSGTTSAT QPSETFMPQRTPSLISGIPTPPSLPGNSKMEAYGGQLGSFPTSQFYVMQAGNA ASSSSSPHMFGGSHMQSSYNAFSLHNPYNLYGYNFPTSPRLAASPEKLSASQS TLLCSSPSNGAFGERQYLPSGMEHSMHMISPSPNNQATNTCDGRQYGAVPGS SSQMSVH
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	70.4 kDa
Uniprot ID	Q965F7
Expiration Date	6 months from date of receipt.