

Product Datasheet

Human LMAN1L protein (orb244943)

Description	Recombinant human Protein ERGIC-53-like
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	GCPPLRRFEYKLSFKGPRLALPGAGIPFWSHHGDAILGLEEVRLTPSMRNRSGAV WSRASVPFSAWEVEVQMRVTGLGRRGAQGMVWYTRGRGHVGSVLGGLAS WDGIGIFFDSPAEDTQDSPAIRVLASDGHIPSEQPGDGASQGLGSCHWDFRNR HPFRARITYWQQLRMSLNGLTPSDPGEFCVDVGPLLLVPGGFFGVSAATGTL ADDHDVLSFLTFSLSEPSPEVPPQPFLEMQQLRLARQLEGLWARLGLGTREDVT PKSDSEAQGEGERLFDLEETLGRHRRILQALRGLSKQLAQERQWKKQLGPPG QARPDGGWALDASCQIPSTPGRGGHLSMSLNKDSAKVGALLHGQWTLQLALQ EMRDAAVRMAAEQVSYLPVGIEHHFLELDHILGLLQEELRGPAAAAKAPRPPG QPPRASSCLQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	63.5 kDa
Uniprot ID	Q9HAT1
Expiration Date	6 months from date of receipt.