

Product Datasheet

Human ASH2L protein (orb244978)

Description	Recombinant human Set1/Ash2 histone methyltransferase complex subunit ASH2
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	This is His-SUMO-tag protein
Protein Sequence	MDTQAGSVDEENGRQLGEVELQCGICTKWFTADTFGIDTSSCLPFMTNYSFHC NVCHHSGNTYFLRKQANLKEMCLSANLTWQSRTQDEHPKTMFSKDKDIIPFI DKYWECMTTRQRP GKMTWPNNIVKTM SKERDVFLVKEHPDPGSKDPEEDYPKF GLLDQDLSNIGPAYDNQKQSSAVSTSGNLNGGIAAGSSGKGRGAKRKQQDGGT TGTTKKARSDPLFSAQRLPPHGYPLEHPFNKDG YRYILAEPDPHAPDPEKLELDC WAGKPIPGDLYRACLYERVLLALHDRAPQLKISDDRLTVVGEKGYSMVRASHGV RKGAWYFEITVDEMPPTAARLGWSQPLGNLQAPLG YDKFSYSWRSKKGTFKH QSIGKHYSYG YQGQDVLGFYINLPEDTETAKSLPDTYKDKALIKFKSYLYFEEDF VDKAEKSLKQTPHSEIIFYKNGVNQGVAYKDIFEGVYFPAISLYKSCTVSINFGPCF KYPPKDLTYRPMSDMGWGAVVEHTLADVLYHVETEVDGRRSPPWEP
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	76.2 kDa
Uniprot ID	Q9UBL3
Expiration Date	6 months from date of receipt.