

Product Datasheet

Human MTMR7 protein (orb245000)

Description

Recombinant human Myotubularin-related protein 7

Endotoxins

Not test.

Preservatives

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Form/Appearance

Liquid or Lyophilized powder

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Tag

N-terminal 6xHis-SUMO-tagged

Note

For research use only

Application notes

This is His-SUMO-tag protein

Protein Sequence

MEHIRTPKVENVRLVDRVSPKKAALGTLTLTATHVIFVENSPDARKETWILHSQIS
TIEKQATTATGCPILLIRCKNFQIIQLIPQERDCHDVYISLIRLARPVKYEELYCFSFN
PMLDKEREQGWVLIDLSEETRMGLPNHYWQLSDVNRDYRVCDSTPTLYVP
KSATAHIIVGSSKFRSRRFPVLSYYYKDNHASICRSSQPLSGFSARCLEDEQMLQ
AIRKANPGSDFVYVVDTRPKLNAMANRAAGKGYENEDNYSNIKFFIGIENIHVM
RNSLQKMLEVCELKSPMSDFLWGLENSGWLRIKAIMDAGIFIAKAVSEEGASV
LVHCSDGWDRTAQVCSVASLLLDPHYRTLKGFVLEIKDWISFGHKFNHRYGNL
DGDPEISPVLDQFIECVWQLMEQFPCAFEFNERFLIHQHHIYSCQFGNFCNSQ
KERRELKIQERTYSLWAHLWKNRADYLNPLFRADHSQTQGTLLHPTTPCNFMYK
FWSGMYNRFKGMQPRQSVTDYLMVKEETQQLLEEELEALEERLEKIQKVQLN
CTKVKSKQSEPSKHSGFSTSDNSIANTPQDYSGNMKSFPSRSPSQGDEDSALILT
QDNLKSSDPDLSANSQESGVEDLSCRSPSGGEHAPSEDSGKDRDSDEAVFLT
A

Purity

Greater than 90% as determined by SDS-PAGE.

Source

E.coli

MW

91.8 kDa

Uniprot ID
[Q9Y216](#)
Expiration Date

6 months from date of receipt.