

Product Datasheet

Human PICK1 protein (orb245242)

Description	Recombinant human PRKCA-binding protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	This is GST-tag protein
Protein Sequence	MFADLDYDIEEDKLGIPVPGKVTLQKDAQNLIGISIGGGAQYCPCLYIVQVFDNT PAALDGTVAAGDEITGVNGRSIKGKTKVEVAKMIQEVKGEVTIHYNKLQADPKQ GMSLDIVLKKVKHRLVENMSSGTADALGLSRAILCNDGLVKRLEELERTAELYKG MTEHTKNLLRAFYELSQTHRAFGDVFSVIGVREPQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	48.9 kDa
Uniprot ID	Q9NRD5
Expiration Date	6 months from date of receipt.