

Product Datasheet

Human HLA-DQA1 protein (orb245607)

Description	Recombinant human HLA class II histocompatibility antigen, DQ alpha 1 chain protein
Target	HLA-DQA1
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	EDIVADHVASYGVNLYQSYGPSGQYTHEFDGDEQFYVDLGRKETVWCLPVLRLQ FRFDPPQFALTNIIVLKHNLNSLIKRSNSTAATNEVPEVTVFSKSPVTLGQPNILICL VDNIFPPVNVITWLSNGHSVTEGVSETSFSLKSDHSFFKISYLTLLPSAEESYDCK VEHWGLDKPLLKHWEPEIPAPMSE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	25.4 kDa
Uniprot ID	P01909
Expiration Date	6 months from date of receipt.